

Recombinant Influenza A virus (H17N10) NA/Neuraminidase Protein, C-His

Catalog No: #AP98340

Package Size: #AP98340-1 50ug #AP98340-2 100ug #AP98340-3 1mg

Description

Product Name	Recombinant Influenza A virus (H17N10) NA/Neuraminidase Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	neuraminidase, NA
Accession No.	AFC35430.1
Target Species	Influenza A virus (A/little yellow-shouldered bat/Guatemala/164/2009 (H17N10))
Calculated MW	48.35 kDa
Target Length	Asn35-Ser442
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.