

Recombinant Influenza A virus (H18N11) HA/Hemagglutinin Protein, N-His

Catalog No: #AP98342

Package Size: #AP98342-1 50ug #AP98342-2 100ug #AP98342-3 1mg

Description

Product Name	Recombinant Influenza A virus (H18N11) HA/Hemagglutinin Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	hemagglutinin, HA
Accession No.	AKC43903.1
Target Species	Influenza A virus (A/dark fruit-eating bat/Bolivia/PBV780-781/2011 (H18N11))
Calculated MW	29.09 kDa
Target Length	Gly53-Ile290
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.