

Recombinant Influenza A virus (H7N7) HA/Hemagglutinin Protein, N-His

Catalog No: #AP98379

Package Size: #AP98379-1 50ug #AP98379-2 100ug #AP98379-3 1mg

Description

Product Name	Recombinant Influenza A virus (H7N7) HA/Hemagglutinin Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	P36346
Target Species	Influenza A virus (strain A/Chicken/Weybridge (H7N7) (Influenza A virus (strain A/FPV/Weybridge (H7N7))
Calculated MW	27.69 kDa
Target Length	Lys58-Glu289
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.