

# Recombinant Influenza B virus (Flu B) HA/Hemagglutinin Protein, C-His

Catalog No: #AP98386

Package Size: #AP98386-1 50ug #AP98386-2 100ug #AP98386-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Influenza B virus (Flu B) HA/Hemagglutinin Protein, C-His
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hemagglutinin, Hemagglutinin HA1 chain, Hemagglutinin HA2 chain, HA
Accession No.	QXN15658.1
Target Species	Influenza B virus (B/New Hampshire/01/2021)
Calculated MW	64.77 kDa
Target Length	Ala15-Thr545
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.