

Recombinant Influenza B virus (Flu D) HE Protein, C-His

Catalog No: #AP98388

Package Size: #AP98388-1 50ug #AP98388-2 100ug #AP98388-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Influenza B virus (Flu D) HE Protein, C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hemagglutinin-esterase-fusion glycoprotein, 3.1.1.53, HE
Accession No.	A0A088CKZ5
Target Species	Influenza D virus (D/bovine/Shandong/Y217/2014)
Calculated MW	36.20 kDa
Target Length	Asp80-Tyr400
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.