

# Recombinant PIV5 F/Fusion glycoprotein F0 Protein, N-His

Catalog No: #AP98389

Package Size: #AP98389-1 50ug #AP98389-2 100ug #AP98389-3 1mg

## Description

|                |                                                                                                                                                                       |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant PIV5 F/Fusion glycoprotein F0 Protein, N-His                                                                                                              |
| Host Species   | E. coli                                                                                                                                                               |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                       |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                       |
| Other Names    | Fusion glycoprotein F0, Fusion glycoprotein F2, Fusion glycoprotein F1, F                                                                                             |
| Accession No.  | P04849                                                                                                                                                                |
| Target Species | Parainfluenza virus 5 (strain W3) (PIV5) (Simian virus 5)                                                                                                             |
| Calculated MW  | 41.00 kDa                                                                                                                                                             |
| Target Length  | Ala126-Ser487                                                                                                                                                         |
| Tag Info       | His-tag                                                                                                                                                               |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                            |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.