

Recombinant PIV5 F/Fusion glycoprotein F0 Protein, C-Strep

Catalog No: #AP98390

Package Size: #AP98390-1 50ug #AP98390-2 100ug #AP98390-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant PIV5 F/Fusion glycoprotein F0 Protein, C-Strep
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Fusion glycoprotein F0, Fusion glycoprotein F2, Fusion glycoprotein F1, F
Accession No.	P04849
Target Species	Parainfluenza virus 5 (strain W3) (PIV5) (Simian virus 5)
Calculated MW	56.51 kDa
Target Length	Leu20-Ser487
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.