

Recombinant Mouse MASP1 Protein, N-His

Catalog No: #AP98396

Package Size: #AP98396-1 50ug #AP98396-2 100ug #AP98396-3 1mg

Description

Product Name	Recombinant Mouse MASP1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Mannose-binding lectin-associated serine protease 1, CRARF, PRSS5, CRARF1, Mannose-binding protein-associated serine protease, Ra-reactive factor serine protease p100, Complement factor MASP-3, RaRF, MASP1, Mannan-binding lectin serine protease 1, Serine protease 5, Complement-activating component of Ra-reactive factor, MASP-1,
Accession No.	P98064
Target Species	Mus musculus (Mouse)
Calculated MW	32.15 kDa
Target Length	Val440-Asn704
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.