

Recombinant Mouse RBP4 Protein, C-Fc

Catalog No: #AP98398

Package Size: #AP98398-1 50ug #AP98398-2 100ug #AP98398-3 1mg

Description

Product Name	Recombinant Mouse RBP4 Protein, C-Fc
Host Species	Mammalian cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Rbp4, Retinol-binding protein 4, Plasma retinol-binding protein, PRBP, RBP
Accession No.	Q00724
Target Species	Mus musculus (Mouse)
Calculated MW	49.72 kDa
Target Length	Glu19-Leu201
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.