

Recombinant Human ZMYND8 Protein, N-His

Catalog No: #AP98399

Package Size: #AP98399-1 50ug #AP98399-2 100ug #AP98399-3 1mg

Description

Product Name	Recombinant Human ZMYND8 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Protein kinase C-binding protein 1, PRKCBP1, RACK7, ZMYND8, CTCL-associated antigen se14-3, Cutaneous T-cell lymphoma-associated antigen se14-3, Rack7, KIAA1125, Zinc finger MYND domain-containing protein 8
Accession No.	Q9ULU4
Target Species	Homo sapiens (Human)
Calculated MW	39.49 kDa
Target Length	Val80-Asn400
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.