

Recombinant Human S100A8 Protein, N-Strep

Catalog No: #AP98400

Package Size: #AP98400-1 50ug #AP98400-2 100ug #AP98400-3 1mg

Description

Product Name	Recombinant Human S100A8 Protein, N-Strep
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	S100 calcium-binding protein A8, Urinary stone protein band A, MRP8, Cystic fibrosis antigen, Calprotectin L1L subunit, MRP-8, p8, Calgranulin-A, Migration inhibitory factor-related protein 8, Leukocyte L1 complex light chain, S100A8, Protein S100-A8, CAGA, CFAG
Accession No.	P05109
Target Species	Homo sapiens (Human)
Calculated MW	12.16 kDa
Target Length	Leu2-Glu93
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.