

Recombinant Mouse STAT3 Protein, N-GST

Catalog No: #AP98405

Package Size: #AP98405-1 50ug #AP98405-2 100ug #AP98405-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse STAT3 Protein, N-GST
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Signal transducer and activator of transcription 3, APRF, STAT3, Acute-phase response factor
Accession No.	P42227
Target Species	Mus musculus (Mouse)
Calculated MW	115 kDa
Target Length	Met1-Met770
Tag Info	GST-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.