

Recombinant Human NSUN5 Protein, C-Strep

Catalog No: #AP98406

Package Size: #AP98406-1 50ug #AP98406-2 100ug #AP98406-3 1mg

Description

Product Name	Recombinant Human NSUN5 Protein, C-Strep
Host Species	Insect Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NOL1R, NOL1-related protein, NSUN5, 28S rRNA (cytosine-C(5))-methyltransferase, WBSCR20A, NOL1/NOP2/Sun domain family member 5, WBSCR20, NSUN5A, Williams-Beuren syndrome chromosomal region 20A protein
Accession No.	Q96P11
Target Species	Homo sapiens (Human)
Calculated MW	47.88 kDa
Target Length	Met1-Arg429
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.