

Recombinant European white birch PPlase/Bet v 7 Protein, N-His

Catalog No: #AP98407

Package Size: #AP98407-1 50ug #AP98407-2 100ug #AP98407-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant European white birch PPlase/Bet v 7 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Peptidyl-prolyl cis-trans isomerase, PPlase, ppiase (CyP), Bet v 7
Accession No.	CAC84116.1
Target Species	Betula pendula (European white birch) (Betula verrucosa)
Calculated MW	20.46 kDa
Target Length	Met1-Ser173
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.