

Recombinant SLEV Capsid protein C/Core protein Protein, N-GST & C-His

Catalog No: #AP98419

Package Size: #AP98419-1 50ug #AP98419-2 100ug #AP98419-3 1mg

Description

Product Name	Recombinant SLEV Capsid protein C/Core protein Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Genome polypeptide, Capsid protein C, Core protein
Accession No.	P09732
Target Species	St. louis encephalitis virus (strain MS1-7)
Calculated MW	38.64 kDa
Target Length	Val13-Arg103
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.