

Recombinant JEV Capsid protein C/Core protein Protein, N-GST & C-His

Catalog No: #AP98426

Package Size: #AP98426-1 50ug #AP98426-2 100ug #AP98426-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant JEV Capsid protein C/Core protein Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Structural polyprotein, Capsid protein C, Core protein
Accession No.	P0DOK8
Target Species	Japanese encephalitis virus (strain M28) (JEV)
Calculated MW	37.93 kDa
Target Length	Lys18-Lys104
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.