

Recombinant PDCoV NSP3/PL-PRO Protein, N-His

Catalog No: #AP98434

Package Size: #AP98434-1 50ug #AP98434-2 100ug #AP98434-3 1mg

Description

Product Name	Recombinant PDCoV NSP3/PL-PRO Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Non-structural protein 3, nsp3, PL2-PRO, Papain-like protease, Papain-like proteinase, PL-PRO
Accession No.	UVJ47558.1
Target Species	Porcine deltacoronavirus
Calculated MW	51.08 kDa
Target Length	Leu939-Glu1383
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.