

Recombinant PDCoV NSP7 Protein, N-GST & C-His

Catalog No: #AP98435

Package Size: #AP98435-1 50ug #AP98435-2 100ug #AP98435-3 1mg

Description

Product Name	Recombinant PDCoV NSP7 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Non-structural protein 7, nsp7
Accession No.	UVJ47558.1
Target Species	Porcine deltacoronavirus
Calculated MW	38.86 kDa
Target Length	Gln3092-Gln3188
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.