

Recombinant Human MZF1 Protein, N-His

Catalog No: #AP98443

Package Size: #AP98443-1 50ug #AP98443-2 100ug #AP98443-3 1mg

Description

Product Name	Recombinant Human MZF1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ZNF42, ZSCAN6, MZF, Myeloid zinc finger 1, Zinc finger and SCAN domain-containing protein 6, Zinc finger protein 42, MZF-1, MZF1
Accession No.	P28698
Target Species	Homo sapiens (Human)
Calculated MW	45.92 kDa
Target Length	Arg354-Glu734
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.