

Recombinant Bovine ISG15 Protein, N-His

Catalog No: #AP98444

Package Size: #AP98444-1 50ug #AP98444-2 100ug #AP98444-3 1mg

Description

Product Name	Recombinant Bovine ISG15 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Interferon-induced 17 kDa protein, G1P2, hUCRP, Ubiquitin-like protein ISG15, IP17, UCRP, ISG15, Ubiquitin cross-reactive protein, Interferon-induced 15 kDa protein,
Accession No.	O02741
Target Species	Bos taurus (Bovine)
Calculated MW	19.47 kDa
Target Length	Met1-Gly154
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.