

Recombinant Human TMED3 Protein, N-His

Catalog No: #AP98449

Package Size: #AP98449-1 50ug #AP98449-2 100ug #AP98449-3 1mg

Description

Product Name	Recombinant Human TMED3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	C15orf22, Membrane protein p24B, p24gamma4, p26, Transmembrane emp24 domain-containing protein 3, p24 family protein gamma-4, TMED3
Accession No.	Q9Y3Q3
Target Species	Homo sapiens (Human)
Calculated MW	13.75 kDa
Target Length	Gly26-Glu123
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.