

Recombinant Human PSME3 Protein, N-His

Catalog No: #AP98450

Package Size: #AP98450-1 50ug #AP98450-2 100ug #AP98450-3 1mg

Description

Product Name	Recombinant Human PSME3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Activator of multicatalytic protease subunit 3, Proteasome activator complex subunit 3, PA28g, 11S regulator complex subunit gamma, PA28gamma, Ki nuclear autoantigen, REG-gamma, PSME3, Proteasome activator 28 subunit gamma
Accession No.	P61289
Target Species	Homo sapiens (Human)
Calculated MW	19.25 kDa
Target Length	Ser111-Tyr254
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.