

Recombinant HCoV-OC43 M/Membrane protein Protein, N-His

Catalog No: #AP98452

Package Size: #AP98452-1 50ug #AP98452-2 100ug #AP98452-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HCoV-OC43 M/Membrane protein Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Membrane protein, M protein, E1 glycoprotein, Matrix glycoprotein, Membrane glycoprotein, M
Accession No.	Q01455
Target Species	Human coronavirus OC43 (HCoV-OC43)
Calculated MW	16.82 kDa
Target Length	Asn103-Ile230 (Val201Ile)
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.