

Recombinant Human TRAM1 Protein, N-His-KSI

Catalog No: #AP98455

Package Size: #AP98455-1 50ug #AP98455-2 100ug #AP98455-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TRAM1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	TRAM, TRAM1, Translocating chain-associated membrane protein 1
Accession No.	Q15629
Target Species	Homo sapiens (Human)
Calculated MW	22.09 kDa
Target Length	Phe318-Ser374
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.