

Recombinant Human CD304/NRP1 Protein, N-His

Catalog No: #AP98457

Package Size: #AP98457-1 50ug #AP98457-2 100ug #AP98457-3 1mg

Description

Product Name	Recombinant Human CD304/NRP1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NRP1, Neuropilin-1, VEGF165R, NRP, CD304, Vascular endothelial cell growth factor 165 receptor, SARS-like Coronaviruses
Accession No.	O14786
Target Species	Homo sapiens (Human)
Calculated MW	35.23 kDa
Target Length	Lys274-Glu584
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.