

Recombinant Human CD304/NRP1 Protein, C-Fc

Catalog No: #AP98458

Package Size: #AP98458-1 50ug #AP98458-2 100ug #AP98458-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD304/NRP1 Protein, C-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NRP1, Neuropilin-1, VEGF165R, NRP, CD304, Vascular endothelial cell growth factor 165 receptor, SARS-like Coronaviruses
Accession No.	O14786
Target Species	Homo sapiens (Human)
Calculated MW	104.47 kDa
Target Length	Met1-Lys644
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.