

Recombinant Human AXL/UFO Protein, C-His

Catalog No: #AP98463

Package Size: #AP98463-1 50ug #AP98463-2 100ug #AP98463-3 1mg

Description

Product Name	Recombinant Human AXL/UFO Protein, C-His
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	UFO, AXL oncogene, Tyrosine-protein kinase receptor UFO, AXL
Accession No.	P30530
Target Species	Homo sapiens (Human)
Calculated MW	49.62 kDa
Target Length	Ala26-Pro449
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.