

Recombinant Human SCARB1/SR-BI/CLA-1/SRB1 Protein, C-Fc

Catalog No: #AP98464

Package Size: #AP98464-1 50ug #AP98464-2 100ug #AP98464-3 1mg

Description

Product Name	Recombinant Human SCARB1/SR-BI/CLA-1/SRB1 Protein, C-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	CD36, CLA1, CD36 antigen-like 1, CD36L1, SRB1, Collagen type I receptor, thrombospondin receptor-like 1, SCARB1, Scavenger receptor class B member 1, CLA-1, SR-BI, CD36 and LIMPII analogous 1, SARS-like Coronaviruses
Accession No.	Q8WTV0
Target Species	Homo sapiens (Human)
Calculated MW	75.07 kDa
Target Length	Pro33-Tyr443
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.