

Recombinant SARS-CoV-2 S1 Protein, C-Fc

Catalog No: #AP98490

Package Size: #AP98490-1 50ug #AP98490-2 100ug #AP98490-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant SARS-CoV-2 S1 Protein, C-Fc
Host Species	Mammalian Cells
Purification	>80% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Spike protein S1, S1 Protein, Spike glycoprotein, S glycoprotein, E2, Peplomer protein, S
Accession No.	P0DTC2
Target Species	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Calculated MW	105 kDa
Target Length	Gln14-Tyr695
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.