

Recombinant SARS-CoV-2 M Protein, N-His

Catalog No: #AP98499

Package Size: #AP98499-1 50ug #AP98499-2 100ug #AP98499-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant SARS-CoV-2 M Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Membrane protein, M, E1 glycoprotein, Matrix glycoprotein, Membrane glycoprotein, M Protein
Accession No.	P0DTC5
Target Species	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Calculated MW	15.84 kDa
Target Length	Arg101-Gln222
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.