

## Recombinant SARS-CoV-2 NSP8 Protein, N-His

Catalog No: #AP98509

Package Size: #AP98509-1 50ug #AP98509-2 100ug #AP98509-3 1mg

## Description

|                |                                                                                                                                                                           |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant SARS-CoV-2 NSP8 Protein, N-His                                                                                                                                |
| Host Species   | E. coli                                                                                                                                                                   |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                           |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| Other Names    | Non-structural protein 8, nsp8, Replicase polyprotein 1ab, pp1ab, ORF1ab polyprotein                                                                                      |
| Accession No.  | P0DTD1                                                                                                                                                                    |
| Target Species | Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)                                                                                                              |
| Calculated MW  | 24.18 kDa                                                                                                                                                                 |
| Target Length  | Ala3943-Gln4140                                                                                                                                                           |
| Tag Info       | His-tag                                                                                                                                                                   |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.