

Recombinant SARS-CoV-2 NSP13 Protein, N-His

Catalog No: #AP98514

Package Size: #AP98514-1 50ug #AP98514-2 100ug #AP98514-3 1mg

Description

Product Name	Recombinant SARS-CoV-2 NSP13 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Helicase, Hel, Non-structural protein 13, nsp13, Replicase polyprotein 1ab, pp1ab, ORF1ab polyprotein
Accession No.	P0DTD1
Target Species	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Calculated MW	69.15 kDa
Target Length	Ala5325-Gln5925
Tag Info	His-tag
Formulation	Lyophilized from a solution in 50mM Tris-HCl, 0.02% NLS, pH 7.5, 5μM ZnCl ₂ , 2mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.