

Recombinant SARS-CoV-2 ORF8/ns8 Protein, N-His-SUMO

Catalog No: #AP98521

Package Size: #AP98521-1 50ug #AP98521-2 100ug #AP98521-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant SARS-CoV-2 ORF8/ns8 Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ORF8 protein, ORF8, Non-structural protein 8, ns8
Accession No.	P0DTC8
Target Species	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Calculated MW	24.76 kDa
Target Length	Phe16-Ile121
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.