

Recombinant Human CD147/BSG Protein, N-His

Catalog No: #AP98523

Package Size: #AP98523-1 50ug #AP98523-2 100ug #AP98523-3 1mg

Description

Product Name	Recombinant Human CD147/BSG Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Hepatoma-associated antigen, BSG, HAb18G, 5F7, TCSF, OK blood group antigen, Collagenase stimulatory factor, CD147, Leukocyte activation antigen M6, Basigin, Extracellular matrix metalloproteinase inducer, EMMPRIN, Tumor cell-derived collagenase stimulatory factor, SARS-like Coronaviruses
Accession No.	P35613
Target Species	Homo sapiens (Human)
Calculated MW	22.47 kDa
Target Length	Gly140-Ala323
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.