

Recombinant Human NRG1/Neuregulin-1 Protein, C-Fc

Catalog No: #AP98524

Package Size: #AP98524-1 50ug #AP98524-2 100ug #AP98524-3 1mg

Description

Product Name	Recombinant Human NRG1/Neuregulin-1 Protein, C-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Glial growth factor, Neu differentiation factor, SMDF, Breast cancer cell differentiation factor p45, Pro-neuregulin-1, membrane-bound isoform, Acetylcholine receptor-inducing activity, ARIA, GGF, HRGA, Pro-NRG1, HGL, HRG, NDF, NRG1, Sensory and motor neuron-derived factor, Heregulin, SARS-like Coronaviruses
Accession No.	Q02297
Target Species	Homo sapiens (Human)
Calculated MW	53.83 kDa
Target Length	Lys22-Lys246
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.