

# Recombinant SARS-CoV S/Spike glycoprotein Protein, C-His

Catalog No: #AP98536

Package Size: #AP98536-1 50ug #AP98536-2 100ug #AP98536-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                |                                                                                                                                                                           |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant SARS-CoV S/Spike glycoprotein Protein, C-His                                                                                                                  |
| Host Species   | Mammalian Cells                                                                                                                                                           |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                           |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                           |
| Other Names    | Nucleoprotein, Nucleocapsid protein, Protein N, NP, NC                                                                                                                    |
| Accession No.  | P59594                                                                                                                                                                    |
| Target Species | Severe acute respiratory syndrome coronavirus (SARS-CoV)                                                                                                                  |
| Calculated MW  | 137.09 kDa                                                                                                                                                                |
| Target Length  | Met1-Lys1193                                                                                                                                                              |
| Tag Info       | His-tag                                                                                                                                                                   |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                           |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.