

Recombinant SARS-CoV-2 NP/Nucleoprotein (Omicron) Protein, N-His

Catalog No: #AP98537

Package Size: #AP98537-1 50ug #AP98537-2 100ug #AP98537-3 1mg

Description

Product Name	Recombinant SARS-CoV-2 NP/Nucleoprotein (Omicron) Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Nucleoprotein, N, Nucleocapsid protein, NC, Protein N, N Protein, Nucleoprotein (Omicron), SARS-like Virus
Accession No.	UFO69287.1
Target Species	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)
Calculated MW	48.57 kDa
Target Length	Met1-Ala416
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.