

Recombinant Human MUC16/CA125 Protein, C-Strep

Catalog No: #AP98551

Package Size: #AP98551-1 50ug #AP98551-2 100ug #AP98551-3 1mg

Description

Product Name	Recombinant Human MUC16/CA125 Protein, C-Strep
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Mucin-16, MUC-16, Ovarian cancer-related tumor marker CA125, CA-125, Ovarian carcinoma antigen CA125, MUC16
Accession No.	Q8WXI7
Target Species	Human
Calculated MW	112.78 kDa
Target Length	Met13360-Gln14347
Tag Info	Strep-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.