

Recombinant Human COL10A1, N-His

Catalog No: #AP98552

Package Size: #AP98552-1 50ug #AP98552-2 100ug #AP98552-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human COL10A1, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Other Names	COL10A1, Collagen alpha-1(X) chain
Accession No.	Q03692
Target Species	Human
Calculated MW	41.77 kDa
Target Length	Thr547-Met680
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.