

Recombinant HBV C/Capsid protein Protein, C-His

Catalog No: #AP98553

Package Size: #AP98553-1 50ug #AP98553-2 100ug #AP98553-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant HBV C/Capsid protein Protein, C-His
Host Species	E.coli
Purification	>95% by SDS-PAGE.
Applications	SDS-PAGE,WB,ELISA,Immunogen,Bioactivity testing in progress
Other Names	Capsid protein, C, Core antigen, Core protein,HBcAg, p21.5,Hepatitis B Virus (HBV)
Accession No.	P03146
Target Species	Hepatitis B virus
Calculated MW	17.92 kDa
Target Length	Met1-Val149
Tag Info	His-tag
Formulation	0.01M PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.