

Recombinant Human ACE2 Protein, HRP

Catalog No: #AP98554

Package Size: #AP98554-1 50ug #AP98554-2 100ug #AP98554-3 1mg

Description

Product Name	Recombinant Human ACE2 Protein, HRP
Host Species	Mammalian Cells
Applications	ELISA, SDS-PAGE, WB
Other Names	Angiotensin-converting enzyme homolog, ACE2, ACEH, Metalloprotease MPROT15, Angiotensin-converting enzyme-related carboxypeptidase, Angiotensin-converting enzyme 2, ACE-related carboxypeptidase
Accession No.	Q9BYF1
Target Species	Homo sapiens (Human)
Calculated MW	85.2 kDa
Target Length	Ser19-Ser740
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.