

Recombinant Mouse MFGE8 Protein, N-His

Catalog No: #AP98567

Package Size: #AP98567-1 50ug #AP98567-2 100ug #AP98567-3 1mg

Description

Product Name	Recombinant Mouse MFGE8 Protein, N-His
Other Names	Lactadherin; MFG-E8; MFGM; MP47; Mfge8; Milk fat globule-EGF factor 8; SED1; Sperm surface protein SP47
Target Species	Mus musculus (Mouse)
Calculated MW	37.45 kDa
Target Length	Leu152-Cys463
Tag Info	N-His
Formulation	Lyophilized from a solution in PBS pH 7.4, 1 mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

>90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.