

beta Actin Antibody HRP Conjugated

Catalog No: #C00044H

Package Size: #C00044H 100ul

Description

Product Name	beta Actin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	MAP conjugated synthetic peptide aa 1-50 375 derived from human beta Actin
Conjugates	HRP
Target Name	beta-Actin
Other Names	BRWS1; PS1TP5BP1; Actin, cytoplasmic 1; Beta-actin; ACTB
Accession No.	Swiss-Prot#P60709NCBI Gene ID60
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-1:50000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.