

MUC16 Antibody FITC Conjugated

Catalog No: #C00069F

Package Size: #C00069F 100ul

Description

Product Name	MUC16 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 12760-12810 22152 derived from human CA125
Conjugates	FITC
Target Name	CA125
Other Names	CA125; Mucin-16; MUC-16; Ovarian cancer-related tumor marker CA125; CA-125; Ovarian carcinoma antigen CA125; MUC16
Accession No.	Swiss-Prot#Q8WXI7NCBI Gene ID94025
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.