

MUC16 Antibody HRP Conjugated

Catalog No: #C00069H

Package Size: #C00069H 100ul

Description

Product Name	MUC16 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 12760-12810 22152 derived from human CA125
Conjugates	HRP
Target Name	CA125
Other Names	CA125; Mucin-16; MUC-16; Ovarian cancer-related tumor marker CA125; CA-125; Ovarian carcinoma antigen CA125; MUC16
Accession No.	Swiss-Prot#Q8WXI7NCBI Gene ID94025
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.