

TGF beta 3 propeptide Antibody HRP Conjugated

Catalog No: #C00077H

Package Size: #C00077H 100ul

Description

Product Name	TGF beta 3 propeptide Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 40-90 412 derived from human TGF-beta 3 propeptide
Conjugates	HRP
Target Name	TGF beta 3 propeptide
Other Names	ARVD; RNHF; ARVD1; TGF-beta3; Transforming growth factor beta-3; TGF-beta-3; TGFB3
Accession No.	Swiss-Prot#P10600NCBI Gene ID7043
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Involved in embryogenesis and cell differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.