

Calcitonin receptor Antibody

Catalog No: #C00098

Package Size: #C00098-1 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Calcitonin receptor Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF
Species Reactivity	Mouse Rat
Specificity	This antibody detects endogenous levels of Calcitonin receptor.
Immunogen Description	KLH conjugated synthetic peptide aa 25-75/516 derived from rat Calcitonin receptor.
Conjugates	Unconjugated
Target Name	Calcitonin receptor; CT-R; C1A C1B; Calcr
Other Names	Calcitonin receptor; CT-R; C1A C1B; Calcr
Accession No.	Gene Name: CalcrSwiss-Prot#P32214Gene ID:116506
Calculated MW	53kDa
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:100-500

Images



Paraformaldehyde-fixed, paraffin embedded (Rat testis); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (CALCR) Polyclonal Antibody (C00098) at 1:500 overnight at 4°C, followed by a conjugated secondary for 20 minutes and DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.