

GAP43 Antibody FITC Conjugated

Catalog No: #C00112F

Package Size: #C00112F 100ul

Description

Product Name	GAP43 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 1-50 227 derived from mouse GAP43
Conjugates	FITC
Target Name	GAP43
Other Names	B-50; Basp2; GAP-43; Neuromodulin; Axonal membrane protein GAP-43; Calmodulin-binding protein P-57; Growth-associated protein 43; Gap43
Accession No.	Swiss-Prot#P06837NCBI Gene ID14432
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

This protein is associated with nerve growth. It is a major component of the motile growth cones" that form the tips of elongating axons. Plays a role in axonal and dendritic filopodia induction.,,P06837,14432,,<https://www.uniprot.org/uniprot/P06837>,<https://www.ncbi.nlm.nih.gov/gene/?term=14432> 30645,30645,C00113H,,CRP Antibody HRP Conjugated,Rabbit,Polyclonal,,IgG,Purified by Protein A.,IHC-P IHC-F,IHC-P=1:50-200 IHC-F=1:50-200,Hu Ms Rt ,,,,CRP,PTX1; C-reactive protein; CRP,Swiss-Prot#P02741NCBI Gene ID1401,,,,,Secreted,,,,,KLH conjugated synthetic peptide aa 140-190 224 derived from human CRP ,HRP,1mg ml,N A,,,,,,,"10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.