

GFAP Antibody FITC Conjugated

Catalog No: #C00142F

Package Size: #C00142F 100ul

Description

Product Name	GFAP Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-80 432 derived from human GFAP
Conjugates	FITC
Target Name	GFAP
Other Names	Glial fibrillary acidic protein; GFAP
Accession No.	Swiss-Prot#P14136NCBI Gene ID2670
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.