

ACEI Antibody FITC Conjugated

Catalog No: #C00146F

Package Size: #C00146F 100ul

Description

Product Name	ACEI Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide of ACEI (C-PTHIKWGD)
Conjugates	FITC
Target Name	ACEI
Other Names	Angiotensin-converting enzyme inhibitory peptide.
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Angiotensin I converting Enzyme Inhibitor, ACEI also designated Bradykinin-potentiating peptide. The peptide: PTHIKWGD. Processed active peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.